

IRAK4 (Phospho-Thr345/Ser346) Antibody

Catalog No: #12812



Package Size: #12812-1 50ul #12812-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IRAK4 (Phospho-Thr345/Ser346) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	IRAK4 (Phospho-Thr345/Ser346) Antibody detects endogenous levels of IRAK4 only when phosphorylated at Thr345/Ser346
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human IRAK4 around the phosphorylation site of Thr345/Ser346.
Target Name	IRAK4 (Phospho-Thr345/Ser346)
Modification	Phospho
Other Names	IRAK4, IPD1, IRAK-4, NY-REN-64, REN64
Accession No.	Swiss-Prot#: Q9NWZ3NCBI Gene ID: 51135
Uniprot	Q9NWZ3
GeneID	51135;
Target Species	human
Calculated MW	60 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only