

## Tri-Methyl-Histone H3 (K10) Polyclonal Antibody

Catalog No: #HW171

Package Size: #hw171-1 50ul #hw171-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Tri-Methyl-Histone H3 (K10) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Tri-Methyl-Histone H3 (K10)
Immunogen Type	protein
Immunogen Description	Synthetic Tri-Methyl peptide from human protein at AA range: 10
Target Name	Tri-Methyl-Histone H3 (K10)
Accession No.	Swiss-Prot:P68431/Q71DI3/P84243/Q6NXT2NCBI Gene ID:8350
Uniprot	P68431
GeneID	8350;8351;8352;8353;8354;8355;8356;8357;8358;8968;
Calculated MW	15 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1B£B'1500-1B£B'12000

## Images



Western blot analysis of 293T lysate, antibody was diluted at 2000. Secondary antibody: diluted at 1:20000

Note: This product is for in vitro research use only