

FoxO4 (Acetyl Lys189) Polyclonal Antibody

Catalog No: #HW175

Package Size: #HW175-1 50ul #HW175-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FoxO4 (Acetyl Lys189) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous FoxO4 when Acetyl occurs at Lys189
Immunogen Type	protein
Immunogen Description	Synthetic Acetyl peptide from human protein at AA range: 189
Target Name	FoxO4 (Acetyl Lys189)
Other Names	Forkhead box protein O4 (Fork head domain transcription factor AFX1)
Accession No.	Swiss-Prot:P98177NCBI Gene ID:4303
Uniprot	P98177
GeneID	4303;
Calculated MW	55 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western blot analysis of 3T3 mouse-kidney KB K562 HeLa lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only