

HNF-4 $\alpha$  (Acetyl Lys106) Polyclonal Antibody

Catalog No: #HW178

Package Size: #HW178-1 50ul #HW178-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	HNF-4 $\alpha$ (Acetyl Lys106) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous HNF-4 $\alpha$ when Acetyl occurs at Lys106
Immunogen Type	protein
Immunogen Description	Synthetic Acetyl peptide from human protein at AA range: 106
Target Name	HNF-4 $\alpha$ (Acetyl Lys106)
Other Names	Hepatocyte nuclear factor 4-alpha (HNF-4-alpha) (Nuclear receptor subfamily 2 group A member 1) (Transcription factor 14) (TCF-14) (Transcription factor HNF-4)
Accession No.	Swiss-Prot:P41235NCBI Gene ID:3172
Uniprot	P41235
GeneID	3172;
Calculated MW	55 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1B£B1500-1B£B12000

## Images



Western blot analysis of mouse-lung mouse-kidney mouse-liver lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000

---

Note: This product is for in vitro research use only