

PRC1(Phospho-Thr481)Polyclonal Antibody

Catalog No: #12399

Package Size: #12399-1 50ul #12399-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PRC1(Phospho-Thr481)Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous PRC1Phospho occurs at Thr481
Immunogen Type	protein
Immunogen Description	Synthetic Phospho peptide from human protein at AA range: 481
Target Name	PRC1(Phospho-Thr481)
Other Names	Protein regulator of cytokinesis 1
Accession No.	Swiss-Prot:O43663NCBI Gene ID:9055
Uniprot	O43663
GeneID	9055;
Calculated MW	72 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western blot analysis of 293T VEC lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only