

## CLECSF6 Polyclonal Antibody

Catalog No: #46747

Package Size: #46747-1 50ul #46747-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	CLECSF6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	CLECSF6 Polyclonal Antibody detects endogenous levels of CLECSF6
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from CLECSF6 at AA range: 161-210
Target Name	CLECSF6
Other Names	C-type lectin domain family 4 member A (C-type lectin DDB27) (C-type lectin superfamily member 6) (Dendritic cell immunoreceptor) (Lectin-like immunoreceptor)
Accession No.	Swiss-Prot:Q9UMR7NCBI Gene ID:50856
Uniprot	Q9UMR7
GeneID	50856;
Calculated MW	28 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details


Western blotting: 1B£B1500-1B£B12000

Immunohistochemistry: IHC: 1:100-1:300


## Images




Western Blot analysis of HEPG2 cells using CLECSF6 Polyclonal Antibody diluted at 1:500. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-colon-cancer, antibody was diluted at 1:200

---

Note: This product is for in vitro research use only