

## RCC1 Polyclonal Antibody

Catalog No: #46749

Package Size: #46749-1 50ul #46749-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	RCC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous RCC1 protein
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 301-350
Target Name	RCC1
Other Names	RCC1 CHC1
Accession No.	Swiss-Prot:P18754NCBI Gene ID:1104
Uniprot	P18754
GeneID	1104;
Calculated MW	45 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details


Western blotting: 1B£B1500-1B£B110000

Immunohistochemistry: IHC: 1:100-1:300


## Images



Western Blot analysis of HEPG2 cells using Antibody diluted at 800. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-skin, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-thyroid, antibody was diluted at 1:200

---

Note: This product is for in vitro research use only