

## RAB11FIP4 Conjugated Antibody

Catalog No: #C29487



Package Size: #C29487-AF350 100ul #C29487-AF405 100ul #C29487-AF488 100ul  
 #C29487-AF555 100ul #C29487-AF594 100ul #C29487-AF647 100ul  
 #C29487-AF680 100ul #C29487-AF750 100ul #C29487-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RAB11FIP4 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	most applications
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human RAB11FIP4 (NP_116321.2).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	RAB11FIP4; FIP4-Rab11; RAB11-FIP4; rab11 family-interacting protein 4
Accession No.	Swiss-Prot#:Q86YS3NCBI Gene ID:84440
Uniprot	Q86YS3
GeneID	84440;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	72kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

## Background

---

The protein encoded by this gene interacts with RAB11 and is thought to be involved in bringing recycling endosome membranes to the cleavage furrow in late cytokinesis. Hypoxic conditions can lead to an upregulation of the encoded protein and enhance the metastatic potential of hepatocellular carcinoma.

---

Note: This product is for in vitro research use only