

Calcitonin Polyclonal Antibody

Catalog No: #46759

Package Size: #46759-1 50ul #46759-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Calcitonin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous Calcitonin
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 80-141
Target Name	Calcitonin
Other Names	Calcitonin [Cleaved into: Calcitonin; Katalcalcin (Calcitonin carboxyl-terminal peptide) (CCP) (PDN-21)]
Accession No.	Swiss-Prot:P01258NCBI Gene ID:796
Uniprot	P06881
GeneID	796;
Calculated MW	16 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western blot analysis of Jurkat U937 Raji mouse-eye lysate, antibody was diluted at 1000. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only