

ASIC1 Polyclonal Antibody

Catalog No: #46764

Package Size: #46764-1 50ul #46764-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ASIC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous ASIC1
Immunogen Description	Synthetic peptide from human protein at AA range: 220-280
Target Name	ASIC1
Other Names	Acid-sensing ion channel 1 (ASIC1) (Amiloride-sensitive cation channel 2, neuronal) (Brain sodium channel 2) (BNaC2)
Accession No.	Swiss-Prot:P78348NCBI Gene ID:41
Uniprot	P78348
GeneID	41;
Calculated MW	55 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1&pound;&ordm;500-1&pound;&ordm;2000

Images



Western blot analysis of mouse-brain and mouse-spinal-cord, antibody was diluted at 1000. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only