

Flt-4 Polyclonal Antibody

Catalog No: #46767

Package Size: #46767-1 50ul #46767-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Flt-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Flt-4
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 640-700
Target Name	Flt-4
Other Names	Vascular endothelial growth factor receptor 3 (VEGFR-3) (EC 2.7.10.1) (Fms-like tyrosine kinase 4) (FLT-4) (Tyrosine-protein kinase receptor FLT4)
Accession No.	Swiss-Prot:P35916NCBI Gene ID:2324
Uniprot	P35916
GeneID	2324;
Calculated MW	170 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details


Western blotting: 1B£B'1500-1B£B'12000

Immunohistochemistry: IHC 1:50-300

Images



Western blot analysis of K562 3T3 lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-placenta, antibody was diluted at 1:200

Note: This product is for in vitro research use only