

Mac-2BP Polyclonal Antibody

Catalog No: #46771

Package Size: #46771-1 50ul #46771-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Mac-2BP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous Mac-2BP
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 100-160
Target Name	Mac-2BP
Other Names	Galectin-3-binding protein (Basement membrane autoantigen p105) (Lectin galactoside-binding soluble 3-binding protein) (Mac-2-binding protein) (MAC2BP) (Mac-2 BP) (Tumor-associated antigen 90K)
Accession No.	Swiss-Prot:Q08380NCBI Gene ID:3959
Uniprot	Q08380
GeneID	3959;
Calculated MW	65 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B'1500-1B£B'12000

Images



Western blot analysis of human blood lysate, antibody was diluted at 1000. Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only