

IL-19 Polyclonal Antibody

Catalog No: #46780

Package Size: #46780-1 50ul #46780-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	IL-19 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Rt
Specificity	The antibody detects endogenous IL-19
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 51-100
Target Name	IL-19
Other Names	Interleukin-19 (IL-19) (Melanoma differentiation-associated protein-like protein) (NG.1)
Accession No.	Swiss-Prot:Q9UHD0NCBI Gene ID:29949
Uniprot	Q9UHD0
GeneID	29949;
Calculated MW	20 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details


Western blotting: 1B£B1500-1B£B12000

Immunohistochemistry: IHC 1:50-300

Images



Western blot analysis of 293T Hela VEC KB mouse-kidney lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

Note: This product is for in vitro research use only