

## GADD34 Polyclonal Antibody

Catalog No: #46783

Package Size: #46783-1 50ul #46783-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GADD34 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous GADD34
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 430-490
Target Name	GADD34
Other Names	Protein phosphatase 1 regulatory subunit 15A (Growth arrest and DNA damage-inducible protein GADD34) (Myeloid differentiation primary response protein MyD116 homolog)
Accession No.	Swiss-Prot:O75807NCBI Gene ID:23645
Uniprot	O75807
GeneID	23645;
Calculated MW	73 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1B£B1500-1B£B12000

## Images



Western blot analysis of MOUSE-BRAIN lysate, antibody was diluted at 1000. Secondary antibody was diluted at 1:20000

---

Note: This product is for in vitro research use only