

NT-3 Polyclonal Antibody

Catalog No: #46794

Package Size: #46794-1 50ul #46794-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | NT-3 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Applications | IHC |
| Species Reactivity | Hu Ms Rt |
| Specificity | The antibody detects endogenous NT-3 |
| Immunogen Type | protein |
| Immunogen Description | Synthetic peptide from human protein at AA range: 180-230 |
| Target Name | NT-3 |
| Other Names | Neurotrophin-3 (NT-3) (HDNF) (Nerve growth factor 2) (NGF-2) (Neurotrophic factor) |
| Accession No. | Swiss-Prot:P20783NCBI Gene ID:4908 |
| Uniprot | P20783 |
| GeneID | 4908; |
| Concentration | 1 mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C |


Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

Note: This product is for in vitro research use only