

Tenascin-R Polyclonal Antibody

Catalog No: #46798

Package Size: #46798-1 50ul #46798-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Tenascin-R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Tenascin-R
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 1270-1350
Target Name	Tenascin-R
Other Names	Tenascin-R (TN-R) (Janusin) (Restrictin)
Accession No.	Swiss-Prot:Q92752NCBI Gene ID:7143
Uniprot	Q92752
GeneID	7143;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-tonsil, antibody was diluted at 1:200

Note: This product is for in vitro research use only