

# PRKCA/PRKCB/PRKCD/PRKCE/PRKCH/PRKCQ Conjugated Antibody

Catalog No: #C30229

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #C30229-AF350 100ul #C30229-AF405 100ul #C30229-AF488 100ul

#C30229-AF555 100ul #C30229-AF594 100ul #C30229-AF647 100ul

#C30229-AF680 100ul #C30229-AF750 100ul #C30229-Biotin 100ul

## Description

Product Name	PRKCA/PRKCB/PRKCD/PRKCE/PRKCH/PRKCQ Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	most applications
Species Reactivity	Hu
Immunogen Description	A synthetic peptide of human PRKCA/PRKCB/PRKCD/PRKCE/PRKCH/PRKCQ.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Accession No.	Swiss-Prot#:P17252/P05771/Q05655/Q02156/P24723/Q04759NCBI Gene ID:5578/5579/5580/5581/5583/5588
Uniprot	P17252
GeneID	5578;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	Refer to figures
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

---

Note: This product is for in vitro research use only