

CD327 Polyclonal Antibody

Catalog No: #46826

Package Size: #46826-1 50ul #46826-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD327 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous CD327
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 71-120
Target Name	CD327
Other Names	Sialic acid-binding Ig-like lectin 6 (Siglec-6) (CD33 antigen-like 1) (CDw327) (Obesity-binding protein 1) (OB-BP1) (CD antigen CD327)
Accession No.	Swiss-Prot:O43699NCBI Gene ID:946
Uniprot	O43699
GeneID	946;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-spleen, antibody was diluted at 1:200

Note: This product is for in vitro research use only