

NNT-1 Polyclonal Antibody

Catalog No: #46837

Package Size: #46837-1 50ul #46837-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NNT-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous NNT-1
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 171-220
Target Name	NNT-1
Other Names	Cardiotrophin-like cytokine factor 1 (B-cell-stimulating factor 3) (BSF-3) (Novel neurotrophin-1) (NNT-1)
Accession No.	Swiss-Prot:Q9UBD9NCBI Gene ID:23529
Uniprot	Q9UBD9
GeneID	23529;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

Note: This product is for in vitro research use only