

## IL-20 Polyclonal Antibody

Catalog No: #46842

Package Size: #46842-1 50ul #46842-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	IL-20 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous IL-20
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 81-130
Target Name	IL-20
Other Names	Interleukin-20 (IL-20) (Cytokine Zcyto10)
Accession No.	Swiss-Prot:Q9NYY1NCBI Gene ID:50604
Uniprot	Q9NYY1
GeneID	50604;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: IHC 1:50-300

## Images



Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:200

---

Note: This product is for in vitro research use only