

ANKRD2 Conjugated Antibody

Catalog No: #C30592



Package Size: #C30592-AF350 100ul #C30592-AF405 100ul #C30592-AF488 100ul
 #C30592-AF555 100ul #C30592-AF594 100ul #C30592-AF647 100ul
 #C30592-AF680 100ul #C30592-AF750 100ul #C30592-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ANKRD2 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	most applications
Species Reactivity	Ms,Rt
Immunogen Description	Recombinant fusion protein of human ANKRD2 (NP_001123453.1).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ANKRD2; ARPP; ankyrin repeat domain 2
Accession No.	Swiss-Prot#:Q9GZV1NCBI Gene ID:26287
Uniprot	Q9GZV1
GeneID	26287;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	42kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

This gene encodes a protein that belongs to the muscle ankyrin repeat protein (MARP) family. A similar gene in rodents is a component of a muscle stress response pathway and plays a role in the stretch-response associated with slow muscle function. Alternative splicing results in multiple transcript variants encoding different isoforms.

Note: This product is for in vitro research use only