

ZAP-70 Polyclonal Antibody

Catalog No: #46858

Package Size: #46858-1 50ul #46858-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | ZAP-70 Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Applications | IHC |
| Species Reactivity | Hu |
| Specificity | The antibody detects endogenous ZAP-70 |
| Immunogen Type | protein |
| Immunogen Description | Synthetic peptide from human protein at AA range: 391-440 |
| Target Name | ZAP-70 |
| Other Names | Tyrosine-protein kinase ZAP-70 (EC 2.7.10.2) (70 kDa zeta-chain associated protein) (Syk-related tyrosine kinase) |
| Accession No. | Swiss-Prot:P43403NCBI Gene ID:7535 |
| Uniprot | P43403 |
| GeneID | 7535; |
| Concentration | 1 mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C |

Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-lymph, antibody was diluted at 1:200

Note: This product is for in vitro research use only