

# Bob 1 Polyclonal Antibody

Catalog No: #46860

Package Size: #46860-1 50ul #46860-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Bob 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous Bob 1
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 201-250
Target Name	Bob 1
Other Names	POU domain class 2-associating factor 1 (B-cell-specific coactivator OBF-1) (BOB-1) (OCA-B) (OCT-binding factor 1)
Accession No.	Swiss-Prot:Q16633NCBI Gene ID:5450
Uniprot	Q16633
GeneID	5450;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: IHC 1:50-300

## Images



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200

---

Note: This product is for in vitro research use only