

DLEU7 Polyclonal Antibody

Catalog No: #46864

Package Size: #46864-1 50ul #46864-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DLEU7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous DLEU7
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 51-100
Target Name	DLEU7
Other Names	Leukemia-associated protein 7 (Deleted in lymphocytic leukemia 7)
Accession No.	Swiss-Prot:Q6UYE1NCBI Gene ID:220107
Uniprot	Q6UYE1
GeneID	220107;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

Note: This product is for in vitro research use only