

NPFF Conjugated Antibody

Catalog No: #C30794



Package Size: #C30794-AF350 100ul #C30794-AF405 100ul #C30794-AF488 100ul
 #C30794-AF555 100ul #C30794-AF594 100ul #C30794-AF647 100ul
 #C30794-AF680 100ul #C30794-AF750 100ul #C30794-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	NPFF Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	most applications
Species Reactivity	Hu,Rt
Immunogen Description	Recombinant fusion protein of human NPFF (NP_003708.1).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	NPFF; FMRFAL; pro-FMRFamide-related neuropeptide FF
Accession No.	Swiss-Prot#:O15130NCBI Gene ID:8620
Uniprot	O15130
GeneID	8620;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	12kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

This gene encodes a member of the FMRFamide related peptide (FARP) family of neuropeptides. The encoded preproprotein is proteolytically processed to generate multiple amidated peptides. These peptides may play a role in the regulation of heart rate and blood pressure and the modulation of morphine-induced antinociception. Patients with hypertension exhibit decreased expression of the encoded protein. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed.

Note: This product is for in vitro research use only