

GSC Polyclonal Antibody

Catalog No: #46871

Package Size: #46871-1 50ul #46871-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GSC Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous GSC
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 207-256
Target Name	GSC
Other Names	Homeobox protein gooseoid
Accession No.	Swiss-Prot:P56915NCBI Gene ID:145258
Uniprot	P56915
GeneID	145258;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-tonsil, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-tonsil, antibody was diluted at 1:200

Note: This product is for in vitro research use only