

Chordin Polyclonal Antibody

Catalog No: #46873

Package Size: #46873-1 50ul #46873-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Chordin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Chordin
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 905-954
Target Name	Chordin
Other Names	Chordin
Accession No.	Swiss-Prot:Q9H2X0NCBI Gene ID:8646
Uniprot	Q9H2X0
GeneID	8646;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-liver-cancer, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-liver-cancer, antibody was diluted at 1:200

Note: This product is for in vitro research use only