

AVP Receptor V3 Polyclonal Antibody

Catalog No: #46878

Package Size: #46878-1 50ul #46878-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AVP Receptor V3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous AVP Receptor V3
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 271-320
Target Name	AVP Receptor V3
Other Names	Vasopressin V1b receptor (V1bR) (AVPR V1b) (AVPR V3) (Antidiuretic hormone receptor 1b) (Vasopressin V3 receptor)
Accession No.	Swiss-Prot:P47901NCBI Gene ID:553
Uniprot	P47901
GeneID	553;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Immunohistochemistry: IHC 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

Note: This product is for in vitro research use only