

## Collagen I Polyclonal Antibody

Catalog No: #46888

Package Size: #46888-1 50ul #46888-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Collagen I Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	This antibody detects endogenous levels of Collagen I
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human Collagen I. at AA range: 468-542
Target Name	Collagen I
Other Names	Collagen alpha-2(I) chain (Alpha-2 type I collagen)
Accession No.	Swiss-Prot:P08123NCBI Gene ID:1278
Uniprot	P08123
GeneID	1278;
Calculated MW	130 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1B£B1500-1B£B12000

## Images



Western Blot analysis of mouse-brain cells using primary antibody diluted at 1:1000(°C overnight). Secondary antibody B£B1Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25B'£C, 1 hour)

---

Note: This product is for in vitro research use only