

Parafibromin Polyclonal Antibody

Catalog No: #46894

Package Size: #46894-1 50ul #46894-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Parafibromin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	This antibody detects endogenous levels of Parafibromin
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human Parafibromin. at AA range: 51-100
Target Name	Parafibromin
Other Names	Parafibromin (Cell division cycle protein 73 homolog) (Hyperparathyroidism 2 protein)
Accession No.	Swiss-Prot:Q6P1J9NCBI Gene ID:79577
Uniprot	Q6P1J9
GeneID	79577;
Calculated MW	61 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of 1,mouse-liver 2,mouse-brain cells using primary antibody diluted at 1:1000(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25BTF£C, 1 hour)

Note: This product is for in vitro research use only