

Integrin  $\beta$ 7 Polyclonal Antibody

Catalog No: #46901

Package Size: #46901-1 50ul #46901-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Integrin $\beta$ 7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	Integrin $\beta$ 7 Polyclonal Antibody detects endogenous levels of Integrin $\beta$ 7
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human Integrin $\beta$ 7. at AA range: 671-720
Target Name	Integrin $\beta$ 7
Other Names	Integrin beta-7 (Gut homing receptor beta subunit)
Accession No.	Swiss-Prot:P26010NCBI Gene ID:3695
Uniprot	P26010
GeneID	3695;
Calculated MW	130 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1B£B1500-1B£B12000

## Images



Western Blot analysis of mouse-heart cells using primary antibody diluted at 1:1000(°C overnight). Secondary antibody B£B1Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25B'££C, 1 hour)

---

Note: This product is for in vitro research use only