

## COL19A1 Polyclonal Antibody

Catalog No: #46906

Package Size: #46906-1 50ul #46906-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	COL19A1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	COL19A1 Polyclonal Antibody detects endogenous levels of COL19A1
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human COL19A1. at AA range: 161-210
Target Name	COL19A1
Other Names	Collagen alpha-1(XIX) chain (Collagen alpha-1(Y) chain)
Accession No.	Swiss-Prot:Q14993NCBI Gene ID:1310
Uniprot	Q14993
GeneID	1310;
Calculated MW	130 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1B£B'1500-1B£B'12000

## Images



Western Blot analysis of HELA cells using primary antibody diluted at 1:2000(°C overnight). Secondary antibodyB£B'Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25B'£C, 1 hour)

---

Note: This product is for in vitro research use only