

Tropomyosin α Polyclonal Antibody

Catalog No: #46908

Package Size: #46908-1 50ul #46908-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Tropomyosin α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Tropomyosin α Polyclonal Antibody detects endogenous levels of Tropomyosin α
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human Tropomyosin α . at AA range: 101-150
Target Name	Tropomyosin α
Other Names	Tropomyosin alpha-1 chain (Alpha-tropomyosin) (Tropomyosin-1)
Accession No.	Swiss-Prot:P09493NCBI Gene ID:7168
Uniprot	P09493
GeneID	7168;
Calculated MW	38 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B'1500-1B£B'12000

Images



Western Blot analysis of 1,mouse-lung 2,mouse-brain 3,mouse-spleen cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B'1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'1£C, 1 hour)

Note: This product is for in vitro research use only