

## AFG3L2 Polyclonal Antibody

Catalog No: #46915

Package Size: #46915-1 50ul #46915-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	AFG3L2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	AFG3L2 Polyclonal Antibody detects endogenous levels of AFG3L2
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human AFG3L2. at AA range: 744-793
Target Name	AFG3L2
Other Names	AFG3-like protein 2 (EC 3.4.24.-) (Paraplegin-like protein)
Accession No.	Swiss-Prot:Q9Y4W6NCBI Gene ID:10939
Uniprot	Q9Y4W6
GeneID	10939;
Calculated MW	88 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1B£B1500-1B£B12000

## Images



Western Blot analysis of 1,mouse-heart 2,293T cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25B'£C, 1 hour)

---

Note: This product is for in vitro research use only