

BRD9 Antibody

Catalog No: #46354

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	BRD9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total BRD9 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic protein corresponding to residues near the C terminal of human BRD9
Target Name	BRD9
Other Names	PRO9856; LAVS3040
Accession No.	Swiss-Prot:Q9H8M2NCBI Gene ID:65980NCBI Protein:BC041590
Uniprot	Q9H8M2
GeneID	65980;
Calculated MW	67 kDa
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:2000

Immunohistochemistry: 1: 50-300

Images




Gel: 8%SDS-PAGE

lysate: 40 µg, Lane 1-3: HEPG2/HeLa/293T cell lysates,

Primary antibody: 46354/BRD9 Antibody) at dilution 1/550

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,

Exposure time: 5 seconds



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46354(BRD9 Antibody) at dilution 1/80, on the right is treated with fusion protein. (Original magnification: x200)

Background

May play a role in chromatin remodeling and regulation of transcription.

Note: This product is for in vitro research use only