

ZNF35 Conjugated Antibody

Catalog No: #C43490



Package Size: #C43490-AF350 100ul #C43490-AF405 100ul #C43490-AF488 100ul
 #C43490-AF555 100ul #C43490-AF594 100ul #C43490-AF647 100ul
 #C43490-AF680 100ul #C43490-AF750 100ul #C43490-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF35 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Species Reactivity	Hu
Immunogen Description	Synthetic peptide of human ZNF35
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	ZNF35
Other Names	HF10; HF.10; Zfp105
Accession No.	Swiss-Prot#: Q9Y2I8NCBI Gene ID: NP_003411
Uniprot	Q9Y2I8
GeneID	22884;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

May be involved in transcriptional regulation. Involved in cell differentiation and/or proliferation.

Note: This product is for in vitro research use only