

C18orf8 Antibody

Catalog No: #46365

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C18orf8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total C18orf8 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic protein corresponding to residues near the C terminal of human C18orf8
Target Name	C18orf8
Other Names	MIC1; Mic-1; HsT2591
Accession No.	Swiss-Prot:Q96DM3NCBI Gene ID:29919NCBI Protein:BC002950
Uniprot	Q96DM3
GeneID	29919;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1: 30-150

Images



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 46365(C18orf8 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x200)

Background

This gene encodes a colon cancer associated protein.

Note: This product is for in vitro research use only