

ADGRL4 Conjugated Antibody

Catalog No: #C43498



Package Size: #C43498-AF350 100ul #C43498-AF405 100ul #C43498-AF488 100ul
 #C43498-AF555 100ul #C43498-AF594 100ul #C43498-AF647 100ul
 #C43498-AF680 100ul #C43498-AF750 100ul #C43498-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	ADGRL4 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Species Reactivity	Hu, Ms, Rt
Immunogen Description	Synthetic peptide of human ADGRL4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	ADGRL4
Other Names	ETL; ELTD1; KPG_003
Accession No.	Swiss-Prot#: Q96RE9NCBI Gene ID: NP_071442
Uniprot	Q96RE9
GeneID	91975;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at -20°C/1 year

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Endothelial orphan receptor that acts as a key regulator of angiogenesis.

Note: This product is for in vitro research use only