

# Rat HB-beta1 ELISA Kit

Catalog No: #EK20108



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat HB-beta1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rt

## Application Details

Detect Range: 0.1 - 16 ng/mL Sample Type: Serum, plasma, Urine, tissue homogenates and other biological fluids.

Sample Volume: 50ul

Assay Time: 3-5h

Detection wavelength: 450nm

Note: This product is for in vitro research use only