

Human Melanocortin receptor 4 (MC4R)ELISA Kit

Catalog No: #EK20115

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Melanocortin receptor 4 (MC4R)ELISA Kit
Brief Description	Elisa Kit
Applications	ELISA
Specificity	Homo sapiens human
Other Names	Melanocortin receptor 4,MC4-R,MC4R
Accession No.	UniProtKB:P32245 GenelDoO 4160
Uniprot	P32245
GeneID	4160;
Storage	12 months

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.060ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 100 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only