

Recombinant human Muscarinic acetylcholine receptor M3

Catalog No: #AP73326

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP73326-1 20ug #AP73326-2 100ug #AP73326-3 1mg

Description

Product Name	Recombinant human Muscarinic acetylcholine receptor M3
Brief Description	Recombinant human Muscarinic acetylcholine receptor M3
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Target Name	CHRM3
Accession No.	Swiss-Prot#:P20309
Uniprot	P20309
GeneID	1131;
Calculated MW	40.7 kDa
Tag Info	N-terminal 6xHis-B2M-tagged
Target Sequence	RIYKETEKRTKELAGLQASGTEAETENFVHPTGSSRSCSSYELQQQSMKRSNRRKYGRCHFWFFTTKSWKPS SEQMDQDHSSSDSWNNNDAAASLENSASSDEEDIGSETRAIYSIVLKLPGHSTILNSTKLPSSDNLQVPEEELG MVDLERKADKLQAQKSVDDGGSPKSFSLPIQLESVDTAKTSDVNSSVGKSTATLPLSFKEATLAKRFALKT RSQITKRKRMSLVKEKAAQT
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Application Details

Tissue Specificity:Homo sapiens (Human)

Expression Region:253-492aa

Sequence Info:Partial

Note: This product is for in vitro research use only