

## CCT8L2 Antibody

Catalog No: #46431

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	CCT8L2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CCT8L2 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic protein corresponding to residues near the C terminal of human CCT8L2
Target Name	CCT8L2
Other Names	CESK1
Accession No.	Swiss-Prot:Q96SF2NCBI Gene ID:150160NCBI Protein:BC033797
Uniprot	Q96SF2
GeneID	150160;
Concentration	1.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1: 30-150

## Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46431(CCT8L2 Antibody) at dilution 1/45, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 46431 (CCT8L2 Antibody) at dilution 1/45, on the right is treated with fusion protein. (Original magnification: x200)

## Background

---

Possible molecular chaperone; assists the folding of proteins upon ATP hydrolysis.

---

Note: This product is for in vitro research use only