

CDK2AP2 Antibody

Catalog No: #46453

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CDK2AP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CDK2AP2 protein.
Immunogen Type	peptide
Immunogen Description	Full length fusion protein of human CDK2AP2
Target Name	CDK2AP2
Other Names	p14; DOC-1R
Accession No.	Swiss-Prot:O75956NCBI Gene ID:10263NCBI Protein:BC002850
Uniprot	O75956
GeneID	10263;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1: 25-100

Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46453(CDK2AP2 Antibody) at dilution 1/20, on the right is treated with fusion protein. (Original magnification: x200)

Background

This gene encodes a protein that interacts with cyclin-dependent kinase 2 associated protein 1. Pseudogenes associated with this gene are located on chromosomes 7 and 9. Alternatively spliced transcript variants have been observed for this gene.?

Note: This product is for in vitro research use only