

CEP89 Antibody

Catalog No: #46479

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CEP89 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CEP89 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide corresponding to internal residues of human CEP89
Target Name	CEP89
Other Names	CEP123; CCDC123
Accession No.	Swiss-Prot:Q96ST8 NCBI Gene ID:84902NCBI Protein:NP_116205
Uniprot	Q96ST8
GeneID	84902;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1: 10-50

Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46479(CEP89 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: x200)

Background

Required for ciliogenesis. Also plays a role in mitochondrial metabolism where it may modulate complex IV activity.

Note: This product is for in vitro research use only