

Mouse Anti-Monkey IgG (H&L) Biotin

Catalog No: #L35064



Package Size: #L35064-1 100 ul #L35064-2 500 ul #L35064-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Anti-Monkey IgG (H&L) Biotin
Brief Description	Conjugated Secondary Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	E, WB, IH, IC
Species Reactivity	Mk
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Monkey IgG. Do not cross reactivity with human IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Monkey IgG
Conjugates	Biotin
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only