

## Goat Anti-Human IgG Fab FITC

Catalog No: #L35067



Package Size: #L35067-1 100 ul #L35067-2 500 ul #L35067-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat Anti-Human IgG Fab FITC
Brief Description	Conjugated Secondary Antibody
Host Species	Goat
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	IF, FC
Species Reactivity	H
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG Fab. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Human IgG Fab
Conjugates	FITC
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

## Application Details

IF (1/250 - 1/2000), FC (1/250 - 1/2000)

Note: This product is for in vitro research use only