

# Rabbit Anti-Cat IgG (H&L) FITC

Catalog No: #L35076



Package Size: #L35076-1 100 ul #L35076-2 500 ul #L35076-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit Anti-Cat IgG (H&L) FITC
Brief Description	Conjugated Secondary Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	IF, FC
Species Reactivity	Ct
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Cat IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Cat IgG
Conjugates	FITC
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

## Application Details

IF (1/250 - 1/2000), FC (1/250 - 1/2000)

Note: This product is for in vitro research use only