

## Rabbit Anti-Mouse IgG (H&amp;L) FITC

Catalog No: #L35089



Package Size: #L35089-1 100 ul #L35089-2 500 ul #L35089-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | Rabbit Anti-Mouse IgG (H&L) FITC   |
| Brief Description     | Conjugated Secondary Antibody  |
| Host Species          | Rabbit   |
| Clonality             | Polyclonal   |
| Purification          | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.   |
| Applications          | IF, FC   |
| Species Reactivity    | M  |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Mouse IgG. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Mouse IgG  |
| Conjugates            | FITC   |
| Formulation           | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide   |

## Application Details

IF (1/250 - 1/2000), FC (1/250 - 1/2000)

Note: This product is for in vitro research use only