

# Mouse Anti-Monkey IgG (H&L) FITC

Catalog No: #L35096



Package Size: #L35096-1 100 ul #L35096-2 500 ul #L35096-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Anti-Monkey IgG (H&L) FITC
Brief Description	Conjugated Secondary Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.
Applications	IF, FC
Species Reactivity	Mk
Specificity	By immunoelectrophoresis and ELISA this antibody reacts specifically with Monkey IgG. Do not cross reactivity with human IgG. No antibody was detected against non immunoglobulin serum proteins.
Immunogen Description	Monkey IgG
Conjugates	FITC
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide

## Application Details

IF (1/250 - 1/2000), FC (1/250 - 1/2000)

Note: This product is for in vitro research use only