

## Goat Anti-Human IgM mu chain - AF405

Catalog No: #L35146



Package Size: #L35146-1 100 ul #L35146-2 500 ul #L35146-3 1 ml

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                                                          |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Goat Anti-Human IgM mu chain - AF405                                                                                                                     |
| Brief Description     | Conjugated Secondary Antibody                                                                                                                            |
| Host Species          | Goat                                                                                                                                                     |
| Clonality             | Polyclonal                                                                                                                                               |
| Purification          | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.                                        |
| Applications          | E, IH, IF, FC                                                                                                                                            |
| Species Reactivity    | H                                                                                                                                                        |
| Specificity           | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgM. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Human IgM mu chain                                                                                                                                       |
| Conjugates            | AF405                                                                                                                                                    |
| Formulation           | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide                                          |

## Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only